

## Product datasheet for **TA363350**

### Hypk Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse HYPK
Specificity:	<b>Expected reactivity:</b> Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	15 kDa
Gene Name:	huntingtin interacting protein K
Database Link:	<a href="#">NP_080594.2</a> <a href="#">Entrez Gene 67693 Mouse</a> <a href="#">Q9CR41</a>
Background:	Has a chaperone-like activity preventing polyglutamine (polyQ) aggregation of HTT. Protects against HTT polyQ-mediated apoptosis in neuronal cells. Required for optimal NAA10-NAA15 complex-mediated N-terminal acetylation (By similarity).
Synonyms:	C15orf63; HSPC136



[View online »](#)

Product images:



Host: Rabbit  
Target Name: Hypk  
Sample Type: Mouse Thymus Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: HYPK  
Sample Tissue: Mouse Thymus lysates  
Antibody Dilution: 1 $\mu$ g/ml