

Product datasheet for **TA363379**

Cox16 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse COX16
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	8 kDa
Gene Name:	cytochrome c oxidase assembly protein 16
Database Link:	NP_001296737.1 Entrez Gene 66272 Mouse Q9CR63-2
Background:	The function of this protein remains unknown.
Synonyms:	C14orf112; HSPC203



[View online »](#)

Product images:



Host: Rabbit
Target Name: Cox16
Sample Type: Mouse Spleen Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: COX16
Sample Tissue: Mouse Spleen lysates
Antibody Dilution: 1ug/ml