

## Product datasheet for **TA363385**

### Ugt2a2 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle terminal region of mouse UGT2A2
Specificity:	<b>Expected reactivity:</b> Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	60 kDa
Gene Name:	UDP glucuronosyltransferase 2 family, polypeptide A2
Database Link:	<a href="#">NP_001019319.1</a> <a href="#">Entrez Gene 552899 Mouse</a> <a href="#">Q6PDD0</a>
Background:	UDPGT is of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.
Synonyms:	Ugt2a2



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: Ugt2a2  
Sample Type: Mouse Testis Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: UGT2A2  
Sample Tissue: Mouse Testis lysates  
Antibody Dilution: 1 $\mu$ g/ml