

Product datasheet for **TA363450**


Tmem87b Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle terminal region of mouse TMEM87B
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	61 kDa
Gene Name:	transmembrane protein 87B
Database Link:	NP_082524.2 Entrez Gene 72477 Mouse Q8BKU8
Background:	May be involved in retrograde transport from endosomes to the trans-Golgi network (TGN).
Synonyms:	FLJ14681

[View online »](#)

Product images:



Host: Rabbit
Target Name: Tmem87b
Sample Type: Mouse Stomach Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: TMEM87B
Sample Tissue: Mouse Stomach lysates
Antibody Dilution: 1ug/ml