

Product datasheet for **TA363473**

Snrnp27 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse SNRNP27
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	21 kDa
Gene Name:	small nuclear ribonucleoprotein 27 (U4/U6.U5)
Database Link:	NP_079941.1 Entrez Gene 66618 Mouse Q8K194-2
Background:	May play a role in mRNA splicing.
Synonyms:	27K; OTTHUMP00000202210; RY1; U4/U6.U5-27K



[View online »](#)

Product images:



Host: Rabbit
Target Name: Snrnp27
Sample Type: Mouse Heart Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: SNRNP27
Sample Tissue: Mouse Heart lysates
Antibody Dilution: 1 μ g/ml