

Product datasheet for **TA363486**

Sntn Rabbit Polyclonal Antibody

Product data:

| | |
|-------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | WB |
| Reactivity: | Mouse |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | The immunogen is a synthetic peptide directed towards the middle region of mouse SNTN |
| Specificity: | Expected reactivity: Mouse |
| Formulation: | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. |
| Concentration: | lot specific |
| Purification: | Affinity purified |
| Conjugation: | Unconjugated |
| Storage: | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability: | Shelf life: one year from despatch. |
| Predicted Protein Size: | 16 kDa |
| Gene Name: | sentan, cilia apical structure protein |
| Database Link: | NP_808292.1 Entrez Gene 218739 Mouse Q8C9X1 |
| Background: | May be a component of the linker structure that bridges the ciliary membrane and peripheral singlet microtubules. |
| Synonyms: | FLJ44379; S100A1L; S100AL; sentan |



[View online »](#)

Product images:



Host: Rabbit
Target Name: Sntn
Sample Type: Mouse Brain Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: SNTN
Sample Tissue: Mouse Brain lysates
Antibody Dilution: 1ug/ml