

## Product datasheet for **TA363491**

### Gtpbp8 Rabbit Polyclonal Antibody


#### Product data:

|                         |   |
|-------------------------|---|
| Product Type:           | Primary Antibodies  |
| Applications:           | WB  |
| Reactivity:             | Mouse   |
| Host:                   | Rabbit  |
| Clonality:              | Polyclonal  |
| Immunogen:              | The immunogen is a synthetic peptide directed towards the middle region of mouse GTPBP8   |
| Specificity:            | <b>Expected reactivity:</b> Mouse   |
| Formulation:            | Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.                                       |
| Concentration:          | lot specific  |
| Purification:           | Affinity purified   |
| Conjugation:            | Unconjugated  |
| Storage:                | For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Stability:              | Shelf life: one year from despatch.   |
| Predicted Protein Size: | 23 kDa  |
| Gene Name:              | GTP-binding protein 8 (putative)  |
| Database Link:          | <a href="#">NP_001152801.1</a><br><a href="#">Entrez Gene 66067 Mouse</a><br><a href="#">Q9CY28-2</a>                                   |
| Background:             | The function of this protein remains unknown.   |
| Synonyms:               | HSPC135   |



[View online »](#)

Product images:



Host: Rabbit  
Target Name: Gtpbp8  
Sample Type: Mouse Pancreas Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: GTPBP8  
Sample Tissue: Mouse Pancreas lysates  
Antibody Dilution: 1 $\mu$ g/ml