

## Product datasheet for **TA363492**

### **Kbtbd4 Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the C terminal region of mouse KBTBD4
Specificity:	<b>Expected reactivity:</b> Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	60 kDa
Gene Name:	kelch repeat and BTB (POZ) domain containing 4
Database Link:	<a href="#">NP_001298045.1</a> <a href="#">Entrez Gene 67136 Mouse</a> <a href="#">Q8R179</a>
Background:	The function of this protein remains unknown.
Synonyms:	BKLHD4; FLJ10450; HSPC252



[View online »](#)

## Product images:



Host: Rabbit  
Target Name: Kbtbd4  
Sample Type: Mouse Brain Lysate  
Antibody Dilution: 1.0 $\mu$ g/ml

Host: Rabbit  
Target Name: KBTBD4  
Sample Tissue: Mouse Brain lysates  
Antibody Dilution: 1 $\mu$ g/ml