

Product datasheet for TA363547


Prl7a2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse PRL7A2
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	27 kDa
Gene Name:	prolactin family 7, subfamily a, member 2
Database Link:	NP_035298.1 Entrez Gene 19114 Mouse Q9R006
Background:	The function of this protein remains unknown.
Synonyms:	PLP-F; Prlpf



Product images:



Host: Rabbit
Target Name: Prl7a2
Sample Type: Rat Testis Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: PRL7A2
Sample Tissue: Rat Testis lysates
Antibody Dilution: 1ug/ml