

Product datasheet for **TA363550**

Zc2hc1c Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of mouse FAM164C
Specificity:	Expected reactivity: Mouse
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	60 kDa
Gene Name:	zinc finger, C2HC-type containing 1C
Database Link:	NP_766002.2 Entrez Gene 72350 Mouse Q8CCG1
Synonyms:	C14orf140; FAM164C



[View online »](#)

Product images:



Host: Rabbit
Target Name: Fam164c
Sample Type: Mouse Testis Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: FAM164C
Sample Tissue: Mouse Testis lysates
Antibody Dilution: 1 ug/ml