

Product datasheet for **TA363589**

LAX1 Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human LAX1
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	35 kDa
Gene Name:	lymphocyte transmembrane adaptor 1
Database Link:	NP_001129662.1 Entrez Gene 54900 Human Q8IWW1-2
Background:	Negatively regulates TCR (T-cell antigen receptor)-mediated signaling in T-cells and BCR (B-cell antigen receptor)-mediated signaling in B-cells.
Synonyms:	FLJ20340; LAX
Protein Families:	Transmembrane



[View online »](#)

Product images:



Host: Rabbit
Target Name: LAX1
Sample Type: HCT116 Cell Lysate
Antibody Dilution: 1.0µg/ml

Host: Rabbit
Target Name: LAX1
Sample Tissue: Human HCT116 Whole Cell lysates
Antibody Dilution: 1 ug/ml