

Product datasheet for **TA363731**

FKBP1B Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	WB
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	The immunogen is a synthetic peptide directed towards the middle region of human FKBP1B
Specificity:	Expected reactivity: Human
Formulation:	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose. <i>Note that this product is shipped as lyophilized powder to China customers.</i>
Concentration:	lot specific
Purification:	Affinity purified
Conjugation:	Unconjugated
Storage:	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Stability:	Shelf life: one year from despatch.
Predicted Protein Size:	11 kDa
Gene Name:	FK506 binding protein 1B
Database Link:	NP_004107.1 Entrez Gene 2281 Human P68106
Background:	The protein encoded by this gene is a member of the immunophilin protein family, which play a role in immunoregulation and basic cellular processes involving protein folding and trafficking. This encoded protein is a cis-trans prolyl isomerase that binds the immunosuppressants FK506 and rapamycin. It is highly similar to the FK506-binding protein 1A. Its physiological role is thought to be in excitation-contraction coupling in cardiac muscle. There are two alternatively spliced transcript variants of this gene encoding different isoforms.



[View online »](#)

Synonyms: FKBP-1B; FKBP-12.6; FKBP1L; FKBP9; FKBP12.6; h-FKBP-12; OTK4; PKBP1L; PPlase; rotamase

Protein Families: Druggable Genome

Product images:



Host: Rabbit
Target Name: FKBP1B
Sample Type: HT1080 Cell Lysate
Antibody Dilution: 1.0 μ g/ml

Host: Rabbit
Target Name: FKBP1B
Sample Tissue: Human HT1080 Whole Cell lysates
Antibody Dilution: 1 μ g/ml