


## Product datasheet for **TA363849**

### Sirpa Mouse Monoclonal Antibody [Clone ID: OX-41]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	OX-41
Applications:	ELISA, IHC
Reactivity:	Rat
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	signal-regulatory protein alpha
Database Link:	<a href="#">Entrez Gene 25528 Rat P97710</a>
Synonyms:	BIT; CD172A; MFR; MYD-1; MYD1; p84; PTPNS1; SHPS-1; SHPS1; SIRP; SIRP-ALPHA-1; Sirp-alpha-2; Sirp-alpha-3; SIRPalha; SIRPalha2

#### Product images:



TA363849, Clone OX-41, rat spleen, frozen section



[View online »](#)