

## Product datasheet for **TA363898**

### ACTH (POMC) Rabbit Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against ANF (1- 28). Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Mouse, Rabbit, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Ser-Leu-Arg-Arg-Ser-Ser-Cys-Phe-Gly-Gly-Arg-Ile- Asp-Arg-Ile-Gly-Ala-Gln-Ser-Gly-Leu-Gly-Cys-Asn-Ser-Phe-Arg-Tyr-OH, (Disulfide bond) coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffered saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the stock solution by freezing aliquots of the stock solution at -20°C or below.
Gene Name:	proopiomelanocortin
Database Link:	<a href="#">P01189</a>
Background:	Atrial Natriuretic Factor (ANF) is a 28 amino acid peptide hormone secreted mainly by the heart atria in response to atrial stretch. ANF acts on the kidney to increase sodium excretion and glomerular filtration rate (GFR), to antagonize renal vasoconstriction, and to inhibit renin secretion. This antibody was generated by immunization of rabbits with Atrial Natriuretic Factor (1- 28)(ms, rb, rt) coupled to a carrier protein.
Synonyms:	ACTH; adrenocorticotropin; alpha-MSH; beta-endorphin; beta-LPH; beta-MSH; CLIP; Corticotropin-lipotropin; gamma-LPH; gamma-MSH; LPH; met-enkephalin; MSH; NPP; OTTHUMP00000119991; OTTHUMP00000200964; POC; pro-ACTH-endorphin; proopiomelanocortin



[View online »](#)

Product images:



Typical titration curve of ANF (1-28)(ms, rb, rt) in a competitive ELISA with this antibody