

Product datasheet for **TA363905**

Corticotropin Releasing Factor (CRH) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Corticotropin Releasing Factor. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user. Please see www.bma.ch for protocols and general information.
Reactivity:	Human, Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Ser-Glu-Glu-Pro-Pro-Ile-Ser-Leu-Asp-Leu-Thr-Phe- His-Leu-Leu-Arg-Glu-Val-Leu-Glu-Met-Ala-Arg-Ala-Glu-Gln-Leu-Ala-Gln- Gln-Ala-His-Ser-Asn-Arg-Lys-Leu-Met-Glu-Ile-Ile-NH ₂ coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	corticotropin releasing hormone
Database Link:	Entrez Gene 1392 Human P06850
Background:	Corticotropin-Releasing Factor (CRF, also known as Corticoliberin) is a 41 amino acid peptide hormone derived from the 196-amino acid preprohormone. It is secreted in the hypothalamus in response to stress, and stimulates both the synthesis and the secretion of adrenocorticotrophic hormone (ACTH). CRF is also synthesized in peripheral tissues and placenta and it is known to suppress appetite and increase subjective feelings of anxiety. This antibody was generated by immunization of rabbits with Corticotropin-Releasing Factor (CRF) coupled to a carrier protein.



[View online »](#)

Synonyms: CRF