

## Product datasheet for **TA363917**

### GHRH Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against GRF. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Tyr-Ala-Asp-Ala-Ile-Phe-Thr-Asn-Ser-Tyr-Arg-Lys- Val-Leu-Gly-Gln-Leu-Ser-Ala-Arg-Lys-Leu-Leu-Gln-Asp-Ile-Met-Ser-Arg- Gln-Gln-Gly-Glu-Ser-Asn-Gln-Glu-Arg-Gly-Ala-Arg-Ala-Arg-Leu-NH <sub>2</sub> coupled to a carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	growth hormone releasing hormone
Database Link:	<a href="#">Entrez Gene 2691 Human P01286</a>
Background:	Growth hormone-releasing factor (GRF), also known as Somatocrinin, is a 44 amino acid peptide hormone produced in the arcuate nucleus of the hypothalamus. GRF stimulates growth hormone production and release by binding to the GRF receptor (GHRHR) on cells in the anterior pituitary. This antibody was generated by immunization of rabbits with GRF coupled to a carrier protein.
Synonyms:	GHRF; GRF; MGC119781; OTTHUMP00000030915; sermorelin; somatocrinin; somatoliberin; Somatorelin



[View online »](#)