

Product datasheet for **TA363921**

Orexin (HCRT) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Bovine, Canine, Mouse, Porcine, Rat
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	hypocretin (orexin) neuropeptide precursor
Database Link:	Entrez Gene 3060 Human O43612
Synonyms:	Hypocretin; NRCLP1; orexin; OX; PPORX; PPOX



[View online »](#)