

## Product datasheet for **TA363929**

### Orexin (HCRT) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Orexin B. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Arg-Ser-Gly-Pro-Pro-Gly-Leu-Gln-Gly-Arg-Leu-Gln- Arg-Leu-Leu-Gln-Ala-Ser-Gly-Asn-His-Ala-Ala-Gly-Ile-Leu-Thr-Met-NH <sub>2</sub> coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	hypocretin (orexin) neuropeptide precursor
Database Link:	<a href="#">Entrez Gene 3060 Human O43612</a>
Background:	Orexin B is a 28-amino acid neuropeptide hormone involved in the regulation a variety of behavioral and physiological processes in response to motivational stimuli. Orexins (Orexin A and Orexin B) are produced by proteolytic cleavage of prepro- orexin. They were so named for their stimulatory role in feeding behavior, but also play important functions in the organism's energy homeostasis and sleep cycles. They exert their effects through G-protein coupled receptors. This antibody was generated by immunization of rabbits with Orexin B coupled to a carrier protein.
Synonyms:	Hypocretin; NRCLP1; orexin; OX; PPORX; PPOX



[View online »](#)