

## Product datasheet for **TA363930**

### Oxytocin neurophysin 1 (OXT) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Oxytocin. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Cys-Tyr-Ile-Gln-Asn-Cys-Pro-Leu-Gly-NH <sub>2</sub> , (Disulfide bond) coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	oxytocin/neurophysin I prepropeptide
Database Link:	<a href="#">Entrez Gene 5020 Human P01178</a>
Background:	Oxytocin is a neuropeptide hormone produced by the hypothalamus and secreted by the pituitary gland. It is derived by enzymatic cleavage from the peptide precursor encoded by the human OXT gene. Oxytocin is released into the bloodstream as a hormone in response to stretching of the cervix and uterus during labor and it plays an important role in the childbirth. This antibody was generated by immunization of rabbits with Oxytocin coupled to a carrier protein.
Synonyms:	MGC126890; MGC126892; OT; OT-NPI



[View online »](#)