

Product datasheet for **TA363932**

Hcrt Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Orexin B. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Mouse, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Arg-Ser-Gly-Pro-Pro-Gly-Leu-Gln-Gly-Arg-Leu-Gln- Arg-Leu-Leu-Gln-Ala-Ser-Gly-Asn-His-Ala-Ala-Gly-Ile-Leu-Thr-Met-NH ₂ coupled to carrier protein.
Formulation:	This polyclonal antibody is supplied as a lyophilized powder. The powder should be rehydrated with 50ml RIA buffer. Upon reconstitution to 50ml total volume, the solution contains 0.1M sodium phosphate buffer (pH 7.4), 0.05M NaCl, 0.1% BSA, 0.01% NaN ₃ , and 0.1% Triton X-100. Store at 4°C. This should ensure antibody stability for approximately one month.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	hypocretin
Database Link:	Entrez Gene 15171 Mouse O55241
Background:	Orexin B is a 28-amino acid neuropeptide hormone involved in the regulation a variety of behavioral and physiological processes in response to motivational stimuli. Orexins (Orexin A and Orexin B) are produced by proteolytic cleavage of prepro- orexin. They were so named for their stimulatory role in feeding behavior, but also play important functions in the organism's energy homeostasis and sleep cycles. They exert their effects through G-protein coupled receptors. This antibody was generated by immunization of rabbits with Orexin B coupled to a carrier protein.



[View online »](#)

Synonyms: Hypocretin; NRCLP1; orexin; OX; PPORX; PPOX