

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA363948

TAC1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Substance P. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Arg-Pro-Lys-Pro-Gln-Gln-Phe-Phe-Gly-Leu-Met-NH2 coupled to carrier protein.
Formulation:	This polyclonal antibody is supplied as a lyophilized powder. The powder should be rehydrated with 50ml of RIA buffer. Upon reconstitution to 50ml total volume, the solution contains 0.1M sodium phosphate buffer (pH 7.4), 0.05M NaCl, 0.1% BSA, 0.01% NaN3, and 0.1% Triton X-100. Store at 4° - 8°C. This should ensure antibody stability for approximately one month.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	tachykinin precursor 1
Database Link:	<u>Entrez Gene 6863 Human</u> <u>P20366</u>
Background:	Substance P is an undecapeptide resulting from the cleavage of Protachykinin-1 (UniProt P20366) into several biologically active peptides. It is mainly secreted by neurons and is involved in biological processes like nociception and inflammation. Substance P mediates interactions between neurons and immune cells, with nerve- derived substance P modulating immune cell proliferation and cytokine production. This antibody was generated by immunization of rabbits with Substance P coupled to a carrier protein.
Synonyms:	NK2; NKA; NKNA; NPK; OTTHUMP00000204986; PPT; TAC2



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US