


Product datasheet for TA363986

Angiotensinogen (AGT) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	angiotensinogen
Database Link:	Entrez Gene 183 Human P01019
Synonyms:	angiotensinogen; ANHU; FLJ92595; FLJ97926; pre-angiotensinogen; SERPINA8

Product images:



Typical titration curve of Angiotensin I (1-10) in a competitive ELISA with this antibody



[View online >](#)