

Product datasheet for **TA364009**

Enkephalin (PENK) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against BAM-12P. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Try-Gly-Gly-Phe-Met-Arg-Arg-Val-Gly-Arg-Pro-Glu- NH ₂ coupled to carrier protein.
Formulation:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	proenkephalin
Database Link:	Entrez Gene 5179 Human P01210
Background:	Bovine Adrenal Medulla Dodecapeptide (BAM-12P) is a 12-amino acid peptide corresponding to the precursor proenkephalin-A. It is an endogenous opioid peptide originally isolated from bovine adrenal medulla, which mimics the effects of opiate drugs, that are pain perception and stress response by exhibiting κ-opioid receptor selectivity. The BAM-12P was described as precursor for Bovine adrenomedullary Met-Enkephalin. This antibody was generated by immunization of rabbits with BAM-12P coupled to a carrier protein.
Synonyms:	preproenkephalin; proenkephalin



[View online »](#)