

Product datasheet for **TA364019**

Nppb Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	Immunohistochemistry, Immunocytochemistry, and other applications.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	This antibody is supplied as a lyophilized powder. Store at 2° to 8°C upon receipt. For long-term storage, store lyophilized powder at -20°C. Once reconstituted, the materials should be stored frozen in aliquot at -20°C to -60°C. Avoid repeated freeze/thaw cycles.
Gene Name:	natriuretic peptide B
Database Link:	Entrez Gene 25105 Rat P13205
Synonyms:	BNP



[View online »](#)