

Product datasheet for TA364031

Calcitonin (CALCA) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Calcitonin. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Cys-Gly-Asn-Leu-Ser-Thr-Cys-Met-Leu-Gly-Thr-Tyr- Thr-Gln-Asp-Phe-Asn-Lys-Phe-His-Thr-Phe-Pro-Gln-Thr-Ala-Ile-Gly-Val- Gly-Ala-Pro-NH ₂ , (Disulfide bond) coupled to carrier protein.
Formulation:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	calcitonin related polypeptide alpha
Database Link:	Entrez Gene 796 Human P01258
Background:	Calcitonin is a 32 amino acid peptide hormone that defines the calcitonin-like protein family. It is formed by the proteolytic cleavage of a larger prepropeptide, which is the product of the CALC1 gene and secreted by parafollicular cells of the thyroid gland in humans, and in the ultimopharyngeal body in animals. The hormone participates in calcium and phosphorus metabolism by reducing blood calcium, and counteracts the effects of parathyroid hormone and Vitamin D3. This antibody was generated by immunization of rabbits with Calcitonin coupled to a carrier protein.
Synonyms:	CALC1; CGRP; CGRP-I; CGRP1; CT; katalcalcin; KC; MGC126648; PCT; Procalcitonin



[View online »](#)

Product images:



Typical titration curve of Calcitonin in a competitive ELISA with this antibody