

Product datasheet for TA364050

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

VIP Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA

Recommended Dilution: This antibody has been tested and validated in ELISA against Vasoactive Intestinal Peptide.

Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well.

Optimal dilutions should be determined by the end user.

Reactivity: Human, Bovine, Porcine, Rat

Host: Rabbit

Clonality: Polyclonal

Immunogen: Synthetic peptide H-His-Ser-Asp-Ala-Val-Phe-Thr-Asp-Asn-Tyr-Thr-Arg- Leu-Arg-Lys-Gln-Met-

Ala-Val-Lys-Lys-Tyr-Leu-Asn-Ser-Ile-Leu-Asn-NH2 coupled to carrier protein.

Formulation: Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl

distilled water. This will give the equivalent of undiluted antiserum.

Concentration: N/A

Conjugation: Unconjugated

Storage: Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing

and freezing of the antiserum by freezing aliquots at -20°C or below.

Gene Name: vasoactive intestinal peptide

Database Link: Entrez Gene 7432 Human

P01282

Background: Vasoactive intestinal peptide (VIP) is a 28 amino acid neuropeptide that is vasoactive in the

intestine. It belongs to a glucagon/secretin superfamily, the ligand of class II G protein-

coupled receptors. VIP is produced in many tissues of vertebrates including the gut, pancreas,

and suprachiasmatic nuclei of the hypothalamus in the brain. VIP functions as


neuromodulator and neurotransmitter, and it is known to regulate muscle activity, epithelial cell secretion, and blood flow in the gastrointestinal tract. This antibody was generated by immunization of guinea pigs with Vasoactive Intestinal Peptide coupled to a carrier protein.

Synonyms: MGC13587; PHM27





Product images:



Typical titration curve of VIP in a competitive ELISA with the antibody from rabbit host