

## Product datasheet for **TA364058**

### ProDynorphin (PDYN) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	prodynorphin
Database Link:	<a href="#">Entrez Gene 5173 Human P01213</a>
Synonyms:	beta-neoendorphin; beta-neoendorphin-dynorphin; dynorphin; leu-enkephalin; leumorphin; MGC26418; PENKB; preprodynorphin; prodynorphin; rimorphin



[View online »](#)