

Product datasheet for **TA364081**

MLN Rabbit Polyclonal Antibody

Product data:

| | |
|-----------------------|--|
| Product Type: | Primary Antibodies |
| Applications: | ELISA |
| Recommended Dilution: | This antibody has been tested and validated in ELISA against Motilin. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user. |
| Reactivity: | Canine |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic peptide H-Phe-Val-Pro-Ile-Phe-Thr-Tyr-Gly-Glu-Leu-Gln-Arg- Met-Gln-Glu-Lys-Glu-Arg-Asn-Lys-Gly-Gln-NH ₂ coupled to carrier protein. |
| Formulation: | Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum. |
| Concentration: | N/A |
| Conjugation: | Unconjugated |
| Storage: | Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below. |
| Database Link: | P19863 |
| Background: | Motilin is a 22-amino acid polypeptide hormone secreted by endocrine Mo cells in the small intestine, especially duodenum and jejunum. It is known for its ability to stimulate gastric activity; in fact, Motilin is released into the general circulation in humans during the inter-digestive state and is the most important factor in controlling the inter-digestive migrating contractions. This antibody was generated by immunization of rabbits with Motilin coupled to a carrier protein. |
| Synonyms: | MGC138519; motilin |



[View online »](#)