

Product datasheet for **TA364084**

Endothelin 1 (EDN1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Mouse, Rabbit, Rat
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	endothelin 1
Database Link:	Entrez Gene 1906 Human P05305
Synonyms:	ET1; HDLCQ7; PPET1; Preproendothelin-1



[View online »](#)