

Product datasheet for TA364096

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GAL Rabbit Polyclonal Antibody

Product data:

Product Type: Primary Antibodies

Applications: ELISA

Recommended Dilution: This antibody has been tested and validated in ELISA against Galanin. Other applications like

immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions

should be determined by the end user.

Reactivity: Porcine

Host: Rabbit

Clonality: Polyclonal

Immunogen: Synthetic peptide H-Gly-Trp-Thr-Leu-Asn-Ser-Ala-Gly-Tyr-Leu-Leu-Gly- Pro-His-Ala-Ile-Asp-Asn-

His-Arg-Ser-Phe-His-Asp-Lys-Tyr-Gly-Leu-Ala- NH2 coupled to carrier protein.

Formulation: This polyclonal antibody is supplied as a lyophilized powder. The powder should be

rehydrated with 50ml of RIA buffer. Upon reconstitution to 50ml total volume, the solution contains 0.1M sodium phosphate buffer (pH 7.4), 0.05M NaCl, 0.1% BSA, 0.01% NaN3, and 0.1% Triton X-100. Store at 4° - 8°C. This should ensure antibody stability for approximately

one month.

Concentration: N/A

Conjugation: Unconjugated

Storage: Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing

and freezing of the antiserum by freezing aliquots at -20°C or below.

Database Link: P07480

Background: Galanin is a 39 amino acid neuropeptide produced from the cleavage of the 123 amino acid

preprogalanin (encoded by GAL gene). It is an inhibitory, hyperpolarizing neuropeptide and as such inhibits neurotransmitter release. Galanin signaling occurs through three G protein-coupled receptors and it is widely expressed in the brain, spinal cord, and gut of humans as well as other mammals. Galanin has been implicated in biological functions such as

nociception, waking and sleep regulation, cognition, feeding, regulation of mood, and regulation of blood pressure. This antibody was generated by immunization of rabbits with

Galanin coupled to a carrier protein.

Synonyms: GAL1; galanin; GALN; GLNN; GMAP; MGC40167

