

## Product datasheet for **TA364106**

### Gastrin (GAST) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Gastrin-1. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide Pyr-Gly-Pro-Trp-Leu-Glu-Glu-Glu-Glu-Ala-Tyr-Gly- Trp-Met-Asp-Phe-NH <sub>2</sub> coupled to carrier protein.
Formulation:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	gastrin
Database Link:	<a href="#">Entrez Gene 2520 Human P01350</a>




[View online »](#)

**Background:**

Gastrin is a peptide hormone released by G cells in the pyloric antrum of the stomach, duodenum, and the pancreas. It stimulates secretion of gastric acid (HCl) by the parietal cells of the stomach and helps in gastric motility. Gastrin binds to cholecystokinin B receptors to stimulate the release of histamines in enterochromaffin-like cells, and it induces the insertion of K<sup>+</sup>/H<sup>+</sup> ATPase pumps into the apical membrane of parietal cells. The encoded polypeptide is preprogastrin, which is enzymatically cleaved to produce progastrin and then gastrin in various forms such as gastrin-34 ("big gastrin"), gastrin-17 ("little gastrin"), gastrin-14 ("minigastrin"). Gastrin-1, also known as Gastrin-17 or "Little Gastrin", is a 17-amino acid peptide secreted by the G cells in the gastric mucosa. It is one of the major bioactive forms of gastrin found in tissue and plasma. This antibody was generated by immunization of rabbits with Gastrin-1 coupled to a carrier protein.

**Synonyms:**

GAS; gastrin

**Product images:**

Typical titration curve of Gastrin-1 in a competitive ELISA with this antibody