

Product datasheet for **TA364114**

GLP1 (GCG) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Glucagon (1-29). Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human, Rat, Porcine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-His-Ser-Gln-Gly-Thr-Phe-Thr-Ser-Asp-Tyr-Ser-Lys- Tyr-Leu-Asp-Ser-Arg-Arg-Ala-Gln-Asp-Phe-Val-Gln-Trp-Leu-Met-Asn- Thr-OH coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	glucagon
Database Link:	Entrez Gene 2641 Human P01275
Background:	Glucagon is a peptide hormone secreted from the pancreatic Islet of Langerhans alpha-cells. It works to raise the concentration of glucose and fatty acids in the bloodstream, and is considered to be the main catabolic hormone of the body. Glucagon (1-29) acts on hepatic enzymes that regulate glucose production and glycogen synthesis. Excessive amounts of circulating glucagon are implicated in the metabolic dysregulation of type 2 diabetes, since such conditions result in hyperglycemia. This antibody was generated by immunization of rabbits with Glucagon (1-29) coupled to a carrier protein.
Synonyms:	GLP1; GLP2; glucagon; GRPP



[View online »](#)