

Product datasheet for **TA364124**

Ghrh Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Growth Hormone Releasing Factor. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Formulation:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	growth hormone releasing hormone
Database Link:	Entrez Gene 29446 Rat P09916
Background:	Growth hormone-releasing factor (GRF), also known as Somatocrinin, is a 44 amino acid peptide hormone produced in the arcuate nucleus of the hypothalamus. GRF stimulates growth hormone production and release by binding to the GRF receptor (GHRHR) on cells in the anterior pituitary. This antibody was generated by immunization of rabbits with Growth Hormone Releasing Factor coupled to a carrier protein.
Synonyms:	GHRF; GRF; MGC119781; OTTHUMP00000030915; sermorelin; somatocrinin; somatoliberin; Somatorelin



[View online »](#)