

Product datasheet for **TA364151**

ACTH (POMC) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	This polyclonal antibody is supplied as a lyophilized powder. Store at 2°- 8°C upon receipt. For long-term storage, store at -20°C. Avoid repeated freeze/thaw cycles.
Gene Name:	proopiomelanocortin
Database Link:	Entrez Gene 5443 Human P01189
Synonyms:	ACTH; adrenocorticotropin; alpha-MSH; beta-endorphin; beta-LPH; beta-MSH; CLIP; Corticotropin-lipotropin; gamma-LPH; gamma-MSH; LPH; met-enkephalin; MSH; NPP; OTTHUMP00000119991; OTTHUMP00000200964; POC; pro-ACTH-endorphin; proopiomelanocortin



[View online »](#)