

## Product datasheet for **TA364161**

### NPY Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	Immunohistochemistry, Immunocytochemistry, and other applications.
Reactivity:	Bovine, Porcine
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	This antibody is supplied as a lyophilized powder. Store at 2° to 8°C upon receipt. For long-term storage, store lyophilized powder at -20°C. Once reconstituted, the materials should be stored frozen in aliquot at -20°C to -60°C. Avoid repeated freeze/thaw cycles.
Database Link:	<a href="#">Entrez Gene 504216 Bovine Q6RUW3</a>
Synonyms:	OTTHUMP00000201946; OTTHUMP00000201947; PYY4



[View online »](#)