

Product datasheet for **TA364164**

Calcitonin N-Terminal Flanking Peptide Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Calcitonin N-Terminal Flanking Peptide. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Ala-Pro-Phe-Arg-Ser-Ala-Leu-Glu-Ser-Ser-Pro-Ala- Asp-Pro-Ala-Thr-Leu-Ser-Glu-Asp-Glu-Ala-Arg-Leu-Leu-Leu-Ala-Ala- Leu-Val-Gln-Asp-Tyr-Val-Gln-Met-Lys-Ala-Ser-Glu-Leu-Glu-Gln-Glu-Gln- Glu-Arg-Glu-Gly-Ser-Ser-Leu-Asp-Ser-Pro-Arg-Ser-OH coupled to carrier protein.
Formulation:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Background:	Calcitonin N-Terminal Flanking Peptide has been shown to stimulate proliferation of normal and neoplastic human osteoblasts at nanomolar concentrations. It has been demonstrated that intracerebroventricular administration of Calcitonin N-Terminal Flanking Peptide significantly decreased food intake and body weight gain for at least 48 h in conscious, freely moving, and unstressed rats fed ad libitum. This antibody was generated by immunization of rabbits with Calcitonin N-Terminal Flanking Peptide coupled to a carrier protein.



[View online »](#)