

Product datasheet for **TA364177**

CHGA Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Porcine
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	This polyclonal antibody is supplied as a lyophilized powder. Store at 2°- 8°C upon receipt. For long-term storage, store at -20°C. Avoid repeated freeze/thaw cycles.
Database Link:	P04404
Synonyms:	CgA; SP-I



[View online »](#)