

Product datasheet for **TA364179**

CHGA Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Pancreastatin. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Porcine
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Gly-Trp-Pro-Gln-Ala-Pro-Ala-Met-Asp-Gly-Ala-Gly- Lys-Thr-Gly-Ala-Glu-Glu-Ala-Gln-Pro-Pro-Glu-Gly-Lys-Gly-Ala-Arg-Glu- His-Ser-Arg-Gln-Glu-Glu-Glu-Glu-Glu-Thr-Ala-Gly-Ala-Pro-Gln-Gly-Leu- Phe-Arg-Gly-NH ₂ coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Database Link:	P04404
Background:	Pancreastatin is a 49-amino acid peptide found in many cells of the neuroendocrine system and it is produced from the proteolytic cleavage of its precursor, chromogranin A. It is localized to insulin-containing beta cells, glucagon-containing alpha cells, and somatostatin-containing D cells of the pancreatic islet. Pancreastatin inhibits release of somatostatin upon glucose stimulation and it may also control carbohydrate metabolism and hyperglycemia. This antibody was generated by immunization of rabbits with Pancreastatin coupled to a carrier protein.
Synonyms:	CgA; SP-I



[View online »](#)