

## Product datasheet for **TA364184**

### Pth Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against pTH (1- 34). Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Ala-Val-Ser-Glu-Ile-Gln-Leu-Met-His-Asn-Leu-Gly- Lys-His-Leu-Ala-Ser-Val-Glu-Arg-Met-Gln-Trp-Leu-Arg-Lys-Lys-Leu-Gln- Asp-Val-His-Asn-Phe-OH (Disulfide bond) coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	parathyroid hormone
Database Link:	<a href="#">Entrez Gene 24694 Rat P04089</a>
Background:	Parathyroid Hormone (pTH) (1-34), is a peptide fragment of the naturally occurring human parathyroid hormone and is secreted from cells of the parathyroid glands and finds its major target cells in bone and kidney. It is an important regulator of calcium and phosphorus metabolism, playing a role in calcium homeostasis of bone, kidney, breast, and placenta. pTH (1-34) administration suppresses cardiovascular calcification and down-regulates aortic osteogenic programs driven by diabetes and dyslipidemia. This antibody was generated by immunization of rabbits with pTH (1-34) coupled to a carrier protein.
Synonyms:	parathormone; Parathyrin; PTH1



[View online »](#)