

## Product datasheet for **TA364195**

### Parathyroid hormone related protein (PTHrP) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Rat
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	parathyroid hormone-like hormone
Database Link:	<a href="#">Entrez Gene 5744 Human P12272</a>
Synonyms:	HHM; MGC14611; osteostatin; PLP; PTH-rP; PTHR; PTHRP



[View online »](#)