

## **Product datasheet for TA364201**

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **VIP Rabbit Polyclonal Antibody**

**Product data:** 

**Product Type:** Primary Antibodies

Applications: ELISA

**Recommended Dilution:** This antibody has been tested and validated in ELISA against PHI-27. Other applications like

immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions

should be determined by the end user.

Reactivity: Porcine

Host: Rabbit

Clonality: Polyclonal

Immunogen: Synthetic peptide H-His-Ala-Asp-Gly-Val-Phe-Thr-Ser-Asp-Phe-Ser-Arg- Leu-Leu-Gly-Gln-Leu-

Ser-Ala-Lys-Lys-Tyr-Leu-Glu-Ser-Leu-Ile-NH2 coupled to a carrier protein.

**Formulation:** Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute

by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffer

saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.

Concentration: N/A

Conjugation: Unconjugated

**Storage:** Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing

and freezing of the antiserum by freezing aliquots at -20°C or below.

Database Link: P01284

**Background:** Peptide Histidine Isoleucine (PHI), is a peptide which functions as a hormone and contains a

composition of 27 amino acids with histidine on the N-terminus and isoleucine on the C-terminus. It is part of family that plays a vital role in the cell growth rate such as in the intestine as well as in brain. PHI-27 was derived from glucagon family called the pituitary adenylate cyclase-activating polypeptide (PACAP) and it has an amino acid sequence homology to vasoactive intestinal peptide, secretin, glucagon, and other growth hormone releasing factor. This peptide is present within the central nervous system that help regulate food consumption behavior, while at peripheral nervous system this peptide accumulates in

the stomach which controls the digestion of food. This antibody was generated by

immunization of rabbits with PHI-27 coupled to a carrier protein.

