

Product datasheet for **TA364222**

Vasopressin (AVP) Rabbit Polyclonal Antibody


Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against (Arg8)- Vasopressin. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Cys-Tyr-Phe-Gln-Asn-Cys-Pro-Arg-Gly-NH ₂ , (disulfide bond) coupled to carrier protein.
Formulation:	Protein A affinity purified from antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 0.2ml distilled water. This stock solution contains 2mg/ml IgG, phosphate buffered saline pH 7.4 (PBS), and 0.02% (w/v) Thimerosal as a preservative.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the stock solution by freezing aliquots of the stock solution at -20°C or below.
Gene Name:	arginine vasopressin
Database Link:	Entrez Gene 551 Human P01185
Background:	(Arg8)-Vasopressin or argipressin is an antidiuretic peptide hormone that plays a key role in maintaining plasma osmolality. The peptide released from the posterior pituitary gland promotes renal tubular reabsorption of water. In high concentration, it also raises blood pressure by inducing moderate vasoconstriction. Argipressin is involved in the control of circadian rhythm, thermoregulation, and ACTH release in the brain. This antibody was generated by immunization of rabbits with (Arg8)-Vasopressin coupled to a carrier protein.
Synonyms:	ADH; ARVP; AVP-NPII; AVRP; neurohypophyseal; VP



[View online »](#)

Product images:



Typical titration curve of (Arg8)-Vasopressin in a competitive ELISA with this antibody