

Product datasheet for **TA364236**

Normal Rabbit Serum, pretitrated Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	E, R, WB
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Reconstitution Method:	Rehydrate the lyophilized neat serum with 50ml of buffer, this should be the same buffer used in the assays. It is sufficient for 500 tests
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Upon receipt, store lyophilized product at 4°C to -20°C. After reconstitution, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.



[View online »](#)