

Product datasheet for **TA364295**

Nicastrin (NCSTN) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	nicastrin
Database Link:	Entrez Gene 23385 Human Q92542
Synonyms:	APH2; ATAG1874; KIAA0253; nicastrin



[View online »](#)