

Product datasheet for **TA364321**

Amyloid Precursor Protein (APP) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human, Mammalian
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	amyloid beta precursor protein
Database Link:	Entrez Gene 351 Human P05067
Synonyms:	A4; AAA; ABPP; AD1; APPI; CTFgamma; CVAP; OTTHUMP00000096096; PN-II; PN2; PreA4



[View online »](#)