

Product datasheet for **TA364364**

Avp Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Rat
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	arginine vasopressin
Database Link:	Entrez Gene 24221 Rat P01186
Synonyms:	ADH; ARVP; AVP-NPII; AVRP; neurohypophyseal; VP



[View online »](#)