

## Product datasheet for **TA364366**

### **PEDF (SERPINF1) Rabbit Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Concentration:	N/A
Conjugation:	Unconjugated
Gene Name:	serpin family F member 1
Database Link:	<a href="#">Entrez Gene 5176 Human P36955</a>
Synonyms:	EPC-1; PEDF; PIG35; Serpin-F1



[View online »](#)