

Product datasheet for **TA364379**

Calcitonin (CALCA) Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against α -CGRP. Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Guinea Pig
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Ala-Cys-Asp-Thr-Ala-Thr-Cys-Val-Thr-His-Arg-Leu- Ala-Gly-Leu-Leu-Ser-Arg-Ser-Gly-Gly-Val-Val-Lys-Asn-Asn-Phe-Val- Pro-Thr-Asn-Val-Gly-Ser-Lys-Ala-Phe-NH ₂ , (Disulfide bond) coupled to carrier protein.
Formulation:	This polyclonal antibody is supplied as a lyophilized powder. The powder should be rehydrated with 50ml of RIA buffer. Upon reconstitution to 50ml total volume, the solution contains 0.1M sodium phosphate buffer (pH 7.4), 0.05M NaCl, 0.1% BSA, 0.01% NaN ₃ , and 0.1% Triton X-100. Store at 4° - 8°C. This should ensure antibody stability for approximately one month.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	calcitonin related polypeptide alpha
Database Link:	Entrez Gene 796 Human P06881



[View online »](#)

Background: α -Calcitonin Gene-Related Peptide (α -CGRP) is a regulatory neuropeptide of 37 amino acids and is formed by the alternative splicing of the calcitonin / CGRP gene located on chromosome 11. It is widely distributed in both the peripheral and central nervous system, predominantly in cell bodies of the dorsal root ganglion (DRG) and it is known to be a potent vasodilator with inotropic and chronotropic effects. This antibody was generated by immunization of guinea pigs with α -Calcitonin Gene- Related Peptide coupled to a carrier protein.

Synonyms: CALC1; CGRP; CGRP-I; CGRP1; CT; katacalcin; KC; MGC126648; PCT; Procalcitonin