

## Product datasheet for **TA364380**

### ACTH (POMC) Guinea Pig Polyclonal Antibody


#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Recommended Dilution:	This antibody has been tested and validated in ELISA against Adrenocorticotrophic Hormone (1-39). Other applications like immunohistochemistry (IHC), FACS or Western Blot may work as well. Optimal dilutions should be determined by the end user.
Reactivity:	Human
Host:	Guinea Pig
Clonality:	Polyclonal
Immunogen:	Synthetic peptide H-Ser-Tyr-Ser-Met-Glu-His-Phe-Arg-Trp-Gly-Lys-Pro- Val-Gly-Lys-Lys-Arg-Arg-Pro-Val-Lys-Val-Tyr-Pro-Asn-Val-Ala-Glu-Asn- Glu-Ser-Ala-Glu-Ala-Phe-Pro-Leu-Glu-Phe-OH (Disulfide bond) coupled to carrier protein.
Formulation:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50µl distilled water. This will give the equivalent of undiluted antiserum.
Concentration:	N/A
Conjugation:	Unconjugated
Storage:	Original vial: at least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the antiserum by freezing aliquots at -20°C or below.
Gene Name:	proopiomelanocortin
Database Link:	<a href="#">Entrez Gene 5443 Human P01189</a>
Background:	Adrenocorticotrophic Hormone (ACTH) is a 39 amino acid peptide hormone produced and secreted by cells of the anterior pituitary gland. It is an important component in therapeutic and diagnostic procedures. This hormone is central to the hypothalamic-pituitary-adrenal axis, a major neuroendocrine system that controls reactions to stress and regulates many body processes, including digestion, the immune system, mood and emotions, and energy metabolism. This antibody was generated by immunization of rabbits with Adrenocorticotrophic Hormone (1-39) coupled to a carrier protein.

[View online »](#)

**Synonyms:**

ACTH; adrenocorticotropin; alpha-MSH; beta-endorphin; beta-LPH; beta-MSH; CLIP; Corticotropin-lipotropin; gamma-LPH; gamma-MSH; LPH; met-enkephalin; MSH; NPP; OTTHUMP00000119991; OTTHUMP00000200964; POC; pro-ACTH-endorphin; proopiomelanocortin

**Product images:**


Typical titration curve of Adrenocorticotrophic Hormone (1-39) in a competitive ELISA with this antibody