

Product datasheet for **TA364400**

Digoxin Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA
Reactivity:	Digoxin
Host:	Rabbit
Isotype:	IgG
Clonality:	Polyclonal
Immunogen:	KLH conjugated Digoxin
Formulation:	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Concentration:	lot specific
Purification:	Antigen affinity purification
Conjugation:	Unconjugated
Storage:	Store at -20°C.
Stability:	1 year
Background:	Digoxin is a beta blocker drug originally derived from the foxglove plant, <i>Digitalis lanata</i> . Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).



[View online »](#)